



RADIOFREKVENTNI NOŽ




- Radiotalasna hirurška jedinica za malu i srednju hirurgiju
- Precizni rez i koagulacija sa minimalnim oštećenjem okolnog tkiva
- Monopolarni i bipolarni režim rada
- Auto start/stop funkcija u bipolarnom režimu
- Memorisanje podešavanja
- Aktivacija preko drške i/ili nožne pedale
- **Radna frekvencija 4 Mhz**
- Digitalna indikacija izlazne snage
- Maksimalna snaga 150W
- Digitalna regulacija izlaznih parametara
- Mogućnost podešavanja vremenskog intervala aplikacije
- „Self - check“ kontrola
- Kontrola kola neutralne elektrode
- Mogućnost korišćenja jednokratnih neutralnih elektroda
- Podešavanje intenziteta akustičnog indikatora
- MDD 93/42/EC Class IIb
- EN60601

SURTRON® FLASH	160 HF
Reference code	10300.201
Max output power CUT (1)	150 W
Max output power BLEND (1)	100 W
Max output power ENHANCED (1)	100 W
Max output power FORCED COAG	100 W
Max output power SOFT COAG	80 W
Max output power BIPOLAR COAG	60 W
Max output power BIPOLAR ABLATION (1)	60 W
Working frequency (1)	4.0 MHz
RF power pulse emission time	10 - 990 ms
RF power pulse interval	1 - 9
Patient circuit	F
Selectable input voltage	115 - 230 Vac
Mains frequency	50 - 60 Hz
Electrical input power	350 VA
Weight	6 kg
Size WxHxD mm	370x144x319

Isporučuje se sa sledećim priborom:

- Drška sa tasterima za aktivaciju
- Drška za mikrohirurške igle
- Elektrode razne - 10 kom.
- Igle za mikro-hirurgiju - 5 kom.
- Elektroda špatula 7 cm - 3 kom.
- Elektroda igla 7 cm - 3 kom.
- Elektroda kugla 7 cm - 3 kom.
- Kabel za neutralnu elektrodu
- Metalna neutralna elektroda
- Neutralna elektroda jednokratna - 4 kom.
- Nožni prekidač, vodootporni
- Mrežni kabel



Monopolarna konfiguracija



Bipolarna konfiguracija



Za rad u bipolarnom režimu potrebno je posebno poručiti:

- bipolarni adapter
- bipolarni kabel
- bipolarni forceps

